VOCO

Futurabond M+

Instructions for use

0482

Product description:

Futurabond M+ is a dental adhesive for use with methacrylate-based restorative, core build-up or luting materials. It can be applied with phosphoric acid etching (selective enamel etching or total-etch) of the dental hard tissue as well as without (self-etch).

Indications:

- [A1] Adhesive for direct restorations of all classes of cavities using light-curing composite / compomer / Ormocer® materials with a methacrylate base.
 - [A2] Adhesive for core build-ups made of light-curing composite materials.
 - [A3] Adhesive for indirect restorations with light-curing luting composites.
 - [B] Intraoral repair of composite, compomer or Ormocer® restorations, ceramic veneers and all-ceramic restorations without an additional primer.
 - [C] Treatment of hypersensitive tooth necks.
 - [D] Protective varnish for glass ionomer cement restorations.
 - [E1] Sealing of cavities prior to amalgam restorations.
 - [E2] Sealing of cavities and core preparations prior to the temporary luting of indirect restorations.
- In combination with **Futurabond M+ Dual Cure Activator**:
- [F1] Adhesive for direct self-curing or dual-curing composite restorations of all classes of cavities and for core build-ups.
 - [F2] Adhesive for indirect restorations using dual-curing or self-curing composites.
 - [G] Luting of root posts with dual-curing or self-curing luting composites.

Method of use:

Futurabond M+ not only ensures a secure bond to natural dental hard tissue, but also to composites, metals and ceramics.

It is necessary to isolate the operating site from moisture. Any contamination of the cavity or material surface with blood or saliva must be prevented. The use of a rubber dam is recommended.

The subsequent stages of application depend on the indication (see sections [A] - [G]). When using light-curing materials, **Futurabond M+** can be applied without the addition of **Futurabond M+ Dual Cure Activator** (see sections [A] - [E]). The **Futurabond M+ Dual Cure Activator** needs to be added only when working with self-curing or dual-curing materials (see sections [F] - [G]).

[A] Direct or indirect restorations with light-curing restorative, core build-up or luting composites

1. Prepare the cavity and core according to the principles of the adhesive restoration technique, and remove any plaque or tartar in the case of non-carious class V cavities. Clean with water. Dry off excess water with a gentle stream of air. Do not overdry the dentine. The objective is a surface with a silky matt appearance. Protect areas near to the pulp with a suitable cavity lining based on calcium hydroxide.
2. Etching of the dental hard tissue (optional)

In self-etch mode, **Futurabond M+** achieves high levels of adhesive strength on enamel and dentine. Prior etching of the dental hard tissue may optionally be carried out. Only enamel which is **unprepared** must always be etched at the beginning of the bonding process (e.g. in the luting of Maryland bridges).

- Etch the surfaces to be etched using a dental acid-etch agent based on phosphoric acid (see applicable instructions for use).

- Aspirate the acid-etch agent, rinse with water for approx. 20 s. Dry off excess moisture with a gentle stream of air. Do not overdry the dentine. The objective is a surface with a silky matt appearance.

3. Mix one drop of **Futurabond M+** and one drop of **Futurabond M+ Dual Cure Activator** in a mixing palette with a disposable applicator for 3 s.

4. Apply the mixed adhesive evenly on the surfaces of the cavity and rub it in for 20 s with a disposable applicator.

5. Dry off the adhesive layer with dry, oil-free air for at least 5 s in order to remove any solvents.

6. Cure the adhesive layer for 10 s using a commercially available polymerisation device (LED or halogen light with an output of > 500 mW/cm²).

7. Depending on the indication, continue to work with the desired material in accordance with the corresponding instructions for use.

[G] Luting of root posts

1. Thoroughly mix one drop of **Futurabond M+** and one drop of **Futurabond M+ Dual Cure Activator** in a mixing palette with a suitable disposable applicator (e.g. **Endo Tim**) for 3 s.
2. Use the applicator to apply the mixed adhesive into the prepared root canal and to the occlusal surface of the core and rub it in for 20 s.

3. Dry off the adhesive layer with dry, oil-free air for at least 5 s in order to remove any solvents.

4. Cure the adhesive layer for 10 s using a commercially available polymerisation device (LED or halogen light with an output of > 500 mW/cm²).

7. Apply the restorative, core build-up or luting material in accordance with the respective instructions for use.

[B] Intraoral repair of fillings or restorations

(made from composite, glass ceramic, zirconium oxide ceramic, aluminium oxide ceramic, metal (incl. amalgam))

1. Thoroughly roughen the surface of the restoration to which the adhesive is to be applied and sandblast it, preferably with aluminium oxide < 50 µm. In the case of ceramics, remove the varnish 1 mm beyond the margin.

2. Clean thoroughly with water spray and dry with moisture- and oil-free air.

3. Continue with the application of **Futurabond M+** immediately afterwards (see [A] 3. - 6.). The use of an additional primer is not necessary.

4. Carry out the repair using a light-curing composite according to the manufacturer's instructions for use.

[C] Treatment of hypersensitive tooth necks

1. Clean the tooth neck in the usual way (e.g. with polishing paste or using a powder jet).

2. Clean thoroughly with water spray and dry with moisture- and oil-free air. Do not overdry. The objective is a surface with a silky matt appearance.

3. Afterwards, continue with the application of **Futurabond M+** (see [A] 3. - 6.).

4. The oxygen inhibition layer should be removed using a cloth or a foam pellet.

[D] Protective varnish for glass ionomer cement restorations

1. If a matrix was used to facilitate shaping, remove the matrix first.

2. Place one drop of **Futurabond M+** on a mixing palette.

3. Apply **Futurabond M+** to the modelled filling with a disposable applicator.

4. Dry off **Futurabond M+** with dry, oil-free air for at least 5 s in order to remove any solvents.

5. Cure **Futurabond M+** on every surface of the restoration using a commercially available polymerisation device for 10 s (LED or light with an output of > 500 mW/cm²).

6. Once the restoration has completely cured, grind the occlusion and seal the ground areas once more using **Futurabond M+** in order to prevent the restoration from drying out.

7. The oxygen inhibition layer should be removed using a cloth or a foam pellet.

[E] Sealing of cavities prior to amalgam restorations or prior to temporary luting

Note: Do not seal cavities with **Futurabond M+** if temporary luting using a methacrylate-based material is scheduled afterwards.

1. Prepare the cavity according to the further treatment planned. Clean with water. Dry off excess water with a gentle stream of air. Do not overdry the dentine. The objective is a surface with a silky matt appearance. Protect areas near to the pulp with a suitable cavity lining based on calcium hydroxide.

2. Place one drop of **Futurabond M+** on a mixing palette.

3. Apply **Futurabond M+** to the cavity surface to be treated with a disposable applicator and rub it in for 20 s.

4. Dry off **Futurabond M+** with dry, oil-free air for at least 5 s in order to remove any solvents.

5. Cure **Futurabond M+** using a commercially available polymerisation device for 10 s (LED or halogen light with an output of > 500 mW/cm²).

6. Depending on the indication, continue to work with the desired material afterwards in accordance with the corresponding instructions for use.

[F] Direct or indirect restorations with self- or dual-curing restorative, core build-up or luting composites

In combination with **Futurabond M+ Dual Cure Activator**, **Futurabond M+** creates a strong bond to dual-curing and self-curing methacrylate-based composites and luting materials.

[A1] Adhesive for direct restorations of all classes of cavities with lichthärtenden Composite-/Compomer-/Ormocer®-Materialien auf Methacrylatbasis.

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