

**Futurabond® M+**

UNIVERSAL ADHESIVE

# Futurabond® M+

## EVERY DROP IS PROFIT

### The right choice of bonding agent for every situation:

In a single bottle, Futurabond M+ combines all that a dentist expects from an adhesive, and for which usually a number of different preparations would be required. It really can be employed universally and makes for fast and high-quality working. Scientific measurements evidence very good adhesion values for Futurabond M+ on enamel and dentine, thus attesting the highest level of reliability within the broadest range of applications for an adhesive. But even in combination with various other materials, such as zirconium oxide substrates, this universal single-bottle adhesive has already been recognised as having excellent adhesion values. In

cases where special indications such as the intraoral repair of all-ceramic restorations have so far required the use of an additional primer, Futurabond M+ alone now demonstrably achieves equal, or even better, results.

To round off its universal versatility, Futurabond M+ allows each user to choose freely whether or not to precondition the dental hard substance with phosphoric acid. Thanks to the use of a new monomer technology, overetching of the dentine is impossible, and post-operative sensitivity is thus reliably prevented.




\* in combination with Futurabond M+ DCA (dual-cure activator) when using self-curing or dual-curing composites


# Futurabond® M+

FLEXIBLE AND RELIABLE

## Direct restoration



## Repairs of all-ceramic restorations



# Futurabond® M+

## UNIVERSAL ADHESIVE

### Indications

Direct and indirect restorations using any light-curing methacrylate-based restorative, core build-up or luting composite

Intraoral repair of fillings, ceramic veneers and all-ceramic restorations without additional primer

Treatment of hypersensitive tooth necks

Protective varnish for glass ionomer cement restorations

Sealing of cavities prior to amalgam restorations or prior to temporary luting

### In combination with Futurabond® M+ DCA (dual-cure activator):

Direct and indirect restorations using any self-curing or dual-curing methacrylate-based restorative, core build-up or luting composite

Luting of root posts with dual-curing or self-curing luting composites

### Advantages

- Can be applied with or without phosphoric acid etching
- Excellently suitable for direct or indirect restorations
- Secure adhesion to various materials such as metals, zirconium or aluminium oxide, as well as silicate ceramics, without additional primer
- Applied in one layer – only 35 seconds total working time

### Presentation

REF 1515	Bottle 5 ml, accessories
REF 1516	Bottle 3 × 5 ml
REF 1518	Futurabond M+ DCA – bottle 2 ml
REF 2246	Endo Tim, application brushes, to be used in the root canal, 50 pcs.
REF 2247	Single Tim, application brushes, 100 pcs.
REF 2315	Mixing palettes, 20 pcs.



VOCO GmbH  
Anton-Flettner-Straße 1-3  
27472 Cuxhaven  
Germany

Freecall: 00 800 44 444 555  
Fax: +49 (0) 4721-719-140

info@voco.com  
www.voco.com

Available from: